



Vacuum -Cum-Pressure Pump



www.acmastechnologies.com | www.acmas.in

**Vacuum -Cum-Pressure Pump
(Diaphragm type)
(Beltless & oil Free)**



ACMAS Technologies Inc.

CE | An ISO 9001:2000 Company

Sales Office: 312-313, Vardhman Capital Mall, L.S.C. 10, Gulabi Bagh, Delhi-110052, INDIA
Tel: +91-11-23654603 (M) +91-9313971681, 9311039044, 9555351619 Fax: +91-11-23654603
Email: acmastechnologies@gmail.com | obromax.delhi@gmail.com
www.acmastechnologies.com | www.acmas.in

ACMAS Technologies Inc.

CE | An ISO 9001:2000 Company



Vacuum -Cum-Pressure Pump



**Vacuum -Cum-Pressure Pump
(Diaphragm type)
(Beltless & oil Free)**



**Vacuum -Cum-Pressure Pump
(Diaphragm type)
(Beltless & oil Free)**

It is a double stage pump to create vacuum through a diaphragm, where continues operations are required without the fear of sample contaminations through oil vapours where even a slightest Environment contamination through oil mist is forbidden. To work on 220/230 volts A.C.

(a) Capacity 25ltrs. / min. air displacement (DOUBLE HEAD)

(b) Capacity 45ltrs. / min air displacement (DOUBLE HEAD)