



Water Still
(Barnstead Pattern)

Water Still
(Barnstead Pattern)

Water Still (Barnstead Pattern)



ACMAS Technologies Inc.

CE An ISO 9001:2000 Company

Sales Office: 312-313, Vardhman Capital Mall, L.S.C. 10, Gulabi Bagh, Delhi-110052, INDIA
Tel: +91-11-23654603 (M) +91-9313971681, 9311039044, 9555351619 Fax: +91-11-23654603
Email: acmastechnologies@gmail.com | obromax.delhi@gmail.com
www.acmastechnologies.com | www.acmas.in

ACMAS Technologies Inc.

CE An ISO 9001:2000 Company



This water distillation unit is designed to give pyrogen free distilled water. It is more economical as the construction is double wall glass wool insulated.

Construction :

Complete stainless steel thick gauge argon arc welded construction. Boiling chamber is double walled and insulated with glass wool to avoid heat losses. Lid with

condenser is easily removable for cleaning of boiling chamber. Heating element with ejection device as safety. The unit is built in constant water level attachment to maintain water level.

These units are available in following capacity

The output is approximate as it is depended on supply voltage.

| CAPACITY (Out-Put Per hour (Approx) | RATING | ELECTRIC SUPPLY |
|-------------------------------------|--------|-----------------|
| 2 Liters | 2 Kw | 240 V |
| 4 Liters | 3 Kw | 240 V |
| 10 Liters | 6 Kw | 240 V |
| 15 Liters | 10 Kw | 440 V |
| 20 Liters | 14 Kw | 440 V |
| 25 Liters | 18 Kw | 440 V |
| 40 Liters | 27 Kw | 440 V |
| 50 Liters | 36 Kw | 440 V |

